


HOW TO USE SIEMENS URINE MULTISTIX 10 SG

1. Dip the entire stick into a urine specimen cup.
2. Wait 60 seconds before reading the protein result.
3. Line up the stick with the testing chart on the side of the bottle.
4. Read and record the urine protein result on urine flowsheet.
5. Bring your urine record sheet to clinic
6. Call the clinic if with protein in the urine for more than 3 days or if there's swelling.

You can cut this section off so the protein tab is at the top.

Only look at the protein tab (4th tab from the top)



PROTEIN	NEGATIVE	TRACE	mg/dL	30 +	100 ++	300 +++	2000 or more ++++
60 seconds							